

## Adipositas weegschaal type 7703.10 / Adiposity Scale

IJkwaardige klasse III / approved Class III

**Bestelnummer/** Order Number

7703.10.001 - incl. RS 232 compoort / with RS 232 interface

7703.10.002 - zonder RS 232 compoort / without RS 232 interface

IJkwaardig CE klasse III, MDD klasse 1

Approved CE class III, MDD class 1

**Weegvermogen /** Weighing Range

0 - 330 kg

**Schaalverdeling /** Division

0 - 200 kg per 0,1 kg 200 - 330 kg per 0,2 kg



De Soehnle Professional Adipositas schaal: de rustgevende zekerheid van een stabiele positie tijdens het wegen.

The Soehnle Professional Adiposity Scale provides a comfortable sense of security during weighing.

### Uniek groot weegplatform met stevige railing

De Soehnle Professional Adipositas weegschaal is door zijn zeer groot weegoppervlak, zijn lage opstap hoogte en zijn stabiele railing dé beste personenweegschaal voor ziekenhuizen, voedingsadviseurs, adipositas en revalidatie klinieken. De apart verhoogde handgrepen op schouderhoogte bieden ook aan zeer zware mensen oriëntatie en zekerheid bij het opstappen en wegen. Het toestel biedt nauwkeurige ijkwaardige weging, is MMD klasse 1 gecertificeerd en heeft eenvoudige bediening.

De afleesunit is 180° draaibaar: gemakkelijk afleesbaar voor zowel de patient als de gebruiker.

### Eigenschappen

- > Groot weegvlak: 89 x 89 cm, hoogte 8,5 cm
- > Weegvermogen tot 330 kg
- > Railing met handgrepen op 148 cm hoogte
- > 180° draaibare gewicht aflees unit
- > IJkwaardige CE klasse III, MDD klasse 1
- > BMI, hold- en dialyse (speciale tarreer-) functie

### Extremely large platform and support rails

The Soehnle Professional Adiposity Scale is especially suited for dietary consultations, adiposity clinics, hospitals and rehabilitation institutes due to its large and low tread and the stable handrails. Special handholds at eye level provide orientation and safety to very adipose patients. In addition to EU Class III and the MPG certification, the high accuracy, the simple operation and the easy to read rotating indicator are further important factors for weighing operations with adipose, handicapped patients.

### Features

- > Approved to Class III, MPG Class 1
- > Platform size: 890 x 890 mm
- > BMI, hold and dialysis functions
- > Digital indicator rotates by 180°
- > Units in kg or lb

Draaibaar aflees display / Rotating display



Stevige railing met verschillende vasthoud mogelijkheden  
Numerous hand position possibilities